

Final Paper - Myopia Control

Ahmed_Drammeh

5/8/2021

Context

The development of an eye towards **myopia** falls under the category known as emmetropization. This process is guided by visual input, and the mechanisms that coordinate this process are not fully understood. It is assumed that emmetropization occurs via an active mechanism by which defocus drives growth of the eye, and that genetic factors and emmetropization both influence the growth of the eye's axis. It has long been assumed that wearing corrective spectacles might possibly hinder the process of emmetropization in young children, and this assumption has been supported in particular by animal studies. However, undercorrection of myopia in humans has been shown to increase the rate of myopic progression.[@shao2015optical]. There has been some research on causal factors involved in the development of myopia. In particular, **statistics show that prolonged near work correlates with the development of myopia**, but it is still unclear whether a causal relationship exists. There are many hypotheses out there which explain the nature and onset and constant evolution of myopia, and one in particular deals with defocusing of the periphery(using micron scale sized lenses on the periphery of a regular spectacle lens). There are some known solutions to retard the progression of myopia in children, but all of them appear to be quite invasive in nature. **This paper's focus, is on analyzing experimental data generated from a potential technological non-invasive (spectacle lenses with microlenses) solution to retard and or eliminate myopic progression in children.** The dataset is produced via a lens fabricating technology called injection molding.

Introduction

The dataset below describes optical measurements of an Injection molded lenses called **Shark**. This spectacle lens is generated by a process called injection molding lithography, which is a replication technique in the manufacturing industry. Any stamp/substrate with fine micro-details engraved/embossed within the surface can be replicated using said technology. The replicated Shark lenses include various characteristics such as front Diopter (describes the base curvature of the front of the lens), radius, cylinder and refractive index. The goal of this report is to evaluate whether or not different molding conditions used during lens generation, will have an impact on the measured global power for the microlenses. Furthermore, the study includes three different grades of Polycarbonate(PC) resin and two types of CNC steel manufactured insert designs (**D19121901 & D_26012101**), which have different geometric dimensions to alter the optical performance of the final product. The analysis for this paper will focus solely on **D26012101**.




Figure 1: *Spectacle lens with 11 concentric rings around the periphery*

Material candidates:

- MEP grade: *CLV1000* - High viscosity(low flow rate resin)
- MEP grade: *CLV4400* - Low viscosity(high flow rate resin)
- Sabic grade: *OQ1022* - Medium viscosity(medium flow rate resin)
- **30 lids from each grade were molded (x 6 process conditions)**
- **15 lids from each process condition were, *PRE.DRIED* @ 120C, to assess impact of drying on the final results.**

Experimental Approach

After a comprehensive recording of several important molding variables such as *Lid (center)CT thickness*, *Fill-time*, *Injection-Pressure*, *Pack+Hold-time*, subsequent optical measurements were also performed after Injection Molding, and these included two key output variables namely *SR2-Optics* and *Jarvis Clearreader*. SR2 optics provides the **cylinder** of a spectacle lens, which is important in vision correction. These variables alongside many others were recorded on a spreadsheet and analyzed to determine whether certain trends can be observed which may point towards good product quality, or whether inferences can be made to inform future processing strategies.

The analysis will take a close look at the optical measurements from the *Jarvis Clear-reader* and see if the microlenses optical performance have any correlation between the various material candidates and their respective process conditioning. The methodology that will be employed in the analyses portion of this paper will include looking at box-plot, scatter plots and histograms as a descriptive tool to make sense of the data. Variables from the Jarvis will mainly highlight the *Global Power Average* for each microlens spanning Rings 1 -11. and relationships between *drying vs non-drying* will be looked at aswell as changes in processing condition. The *Global Power Average* which is the indicator of the effectiveness of the lens to correct vision, will be compared to the design specification to see if the injection molding process can successfully replicate the microlenses and maintain its optical power profile. The dataset can be cumbersome and tedious, but R tidying & cleaning strategies will be utilized to navigate and capture key parts of the data of interest.

Data Cleaning

The name of the dataset is called **Shark**, which is the products code-name aswell. The dataset was imported into R as a **csv google spreadsheet** obtained from a machine operators run of an injection molding

machine. There are a total of **325 observations** and **38 variables** involved in the raw data set. The data first had to be cleaned since several of the columns had column names with a combination of letters, numbers and symbols, which creates a set of unwanted strings and factors. Using the **stringsAsFactors** option during data import prevents the data-frame from converting character strings as factors. Additionally, the global substitution or **gsub()** function aided in converting the character strings to numeric. Below are the series of codes.

```
Shark<-read.csv("Shark_Encap-PMMA_Over_Mold_V4.csv", stringsAsFactors = FALSE, skip=2, header=TRUE)
dim(Shark)
```

```
## [1] 324 39
```

```
str(Shark)
```

```
## 'data.frame': 324 obs. of 39 variables:
## $ Date : chr "2/6/2021" "" "" "" ...
## $ Sample.ID : chr "A16" "A17" "A18" "A19" ...
## $ CT..mm. : chr "1.16" "1.16" "1.16" "1.16" ...
## $ PRE.DRY : chr "CLV1000-No" "CLV1000-No" "CLV1000-No" "CLV1000-No" ...
## $ Mold : chr "R&D 2-Cav" "R&D 2-Cav" "R&D 2-Cav" "R&D 2-Cav" ...
## $ A.SIDE : int 217 217 217 217 217 217 217 217 217 ...
## $ B.SIDE : int 220 220 220 220 220 220 220 220 220 ...
## $ A.SIDE.1 : int 195 195 195 195 195 195 195 195 195 ...
## $ B.SIDE.1 : int 200 200 200 200 200 200 200 200 200 ...
## $ Melt.Temp..F. : int 460 460 460 460 460 460 460 460 460 ...
## $ Fill.time..s. : num 6.52 6.52 6.52 6.53 6.53 6.56 6.56 6.55 6.55 6.52 ...
## $ Cooling.time : int 325 325 325 325 325 325 325 325 325 325 ...
## $ X1st.... : num 5 5 5 5 5 5 5 5 5 ...
## $ X2nd.... : chr "12" "12" "12" "12" ...
## $ X3rd.... : int 8 8 8 8 8 8 8 8 8 ...
## $ X1st..psi..x.1000 : num 4.2 3.9 3.9 3.9 3.9 4 4 4.1 4.1 3.9 ...
## $ X2nd..psi..x.1000 : num 8.28 8.33 8.33 8.38 8.38 8.38 8.38 8.28 8.28 8.28 ...
## $ Inj.P.Setting.... : int 65 65 65 65 65 65 65 65 65 ...
## $ Pack.Set.... : int 58 58 58 58 58 58 58 58 58 ...
## $ Hold.Set.... : int 55 55 55 55 55 55 55 55 55 ...
## $ Pack.Pressure..psi..x.1000: num 15.1 15 15 15 15 15 15 15.1 15.1 ...
## $ Hold.Pressure..psi..x.1000: num 14.2 14.2 14.2 14.2 14.2 ...
## $ Shot.size..mm. : num 70.8 70.7 70.7 70.6 70.6 ...
## $ Cushion..mm. : num 6.62 6.88 6.88 7.04 7.04 7.04 7.04 7.03 7.03 6.81 ...
## $ Diopter.3.5...0.06 : chr "3.54" NA "3.54" "3.54" ...
## $ Cyl...0.06. : chr "0.02" NA "0.03" "0.02" ...
## $ Insert..Design : chr "D_19121901" "D_19121901" "D_19121901" "D_19121901" ...
## $ R.01 : chr "3.577" NA "3.547" "3.438" ...
## $ R.02 : num 3.54 NA 3.52 3.39 3.47 ...
## $ R.03 : num 3.51 NA 3.5 3.38 3.48 ...
## $ R.04 : chr "3.447" NA "3.428" "3.316" ...
## $ R.05 : num 3.35 NA 3.33 3.22 3.33 ...
## $ R.06 : num 3.23 NA 3.21 3.1 3.22 ...
## $ R.07 : num 3.06 NA 3.05 2.92 3.08 ...
## $ R.08 : num 2.86 NA 2.85 2.73 2.87 ...
## $ R.09 : num 2.58 NA 2.57 2.46 2.59 ...
## $ R.10 : num 2.35 NA 2.34 2.27 2.35 ...
## $ R.11 : num 2 NA 1.98 1.93 1.98 ...
## $ X : int NA NA NA NA NA NA NA NA NA ...
```

```
## Switching from character to numeric to setup data visualization
Shark$R.01<-as.numeric(as.character(gsub(",", "", Shark$R.01)))
Shark$R.04<-as.numeric(as.character(gsub(",", "", Shark$R.04)))
Shark$CT..mm.<-as.numeric(as.character(gsub(",", "", Shark$CT..mm.)))
Shark$Diopter.3.5...0.06<-as.numeric(as.character(gsub(",", "", Shark$Diopter.3.5...0.06)))
Shark$Cyl...0.06.<-as.numeric(as.character(gsub(",", "", Shark$Cyl...0.06.)))
Shark$X<-NULL
## Selecting specific rows of interest involved in the analysis. This includes cylinder (Cyl), Rings 1-
library(tidyverse)
Shark1<-select(Shark,4:5, 26:38)
```


Data Tidying & Visualization

After the initial data cleaning step, it was vital to look at some descriptive statistics of the data and see if trends can be identified. In this data, it was important to shift the data from wide to long format, to be able to compare the variables accurately. Below are some of the visualizations that were used to help guide in the understanding of one/more response variables (global Power average, Cylinder) with respect to input variables (Resin grade, drying... etc). Most of the initial visualizations focused on expected trends and behaviors, as well as trends to help consolidate the data further.


```
library(ggplot2)
library(tidyverse)
library(tidyr)
library(dplyr)
Shark.tidy<-gather(Shark1, key="Ring", value="Global_Power_Avg", -1,-2,-3,-4, na.rm=TRUE)
head(Shark.tidy)
```

##	PRE.DRY	Mold	Cyl...0.06.	Insert..Design	Ring	Global_Power_Avg
## 1	CLV1000-No R&D	2-Cav	0.02	D_19121901	R.01	3.577
## 3	CLV1000-No R&D	2-Cav	0.03	D_19121901	R.01	3.547
## 4	CLV1000-No R&D	2-Cav	0.02	D_19121901	R.01	3.438
## 5	CLV1000-No R&D	2-Cav	0.02	D_19121901	R.01	3.513
## 6	CLV1000-No R&D	2-Cav	0.02	D_19121901	R.01	3.493
## 7	CLV1000-No R&D	2-Cav	0.02	D_19121901	R.01	3.546

```
## Some expected relationships. Insert design D26012101 has a higher power. (Global Power average > D19
Shark.tidy %>%
  mutate(Cyl = Cyl...0.06.) %>%
  ggplot(aes(x=Cyl, y=Global_Power_Avg, color=Insert..Design)) + geom_point() + coord_cartesian(ylim=c(
```



```
Shark.tidy %>%
  mutate(Cyl = Cyl...0.06.) %>%
  ggplot(aes(y =Global_Power_Avg, x=Cyl, fill=Insert..Design)) + geom_boxplot() + ggtitle("Lens Cylinder vs Insert Design")
```




No influence on lens cylinder based on Insert design.


```
Shark.tidy %>%
```

```
  mutate(Cyl = Cyl...0.06.) %>%
```


```
  ggplot(aes(x =Insert..Design, y=Cyl, fill=Insert..Design)) + geom_boxplot() + ggtitle("Lens Cylinder v
```




```
## There Should be no influence between Cylinder vs Global Power average when it comes to Resin grade
Shark.tidy %>%
  mutate(Cyl = Cyl...0.06.) %>%
  ggplot(aes(x=Cyl, y=Global_Power_Avg, color=PRE.DRY)) + geom_point() + coord_cartesian(ylim=c(2,5))
```



```
Shark.tidy %>%
  mutate(Cyl = Cyl...0.06.) %>%
  ggplot(aes(x=Cyl, fill=PRE.DRY)) + geom_histogram(binwidth = 0.01) + ggtitle("Lens Cylinder vs Power")
```





```
Shark.tidy %>%
  mutate(Cyl = Cyl...0.06.) %>%
  ggplot(aes(x=Cyl, y=Global_Power_Avg, fill=PRE.DRY)) + geom_boxplot() + coord_cartesian(ylim=c(2,5))
```



```
## Using filter() to selectively focus on the Insert design (D26012101) we hinted on, in the experiment.
Shark_tidy1<-filter(Shark.tidy, Insert..Design=="D_26012101")%>% mutate(Cyl...0.06.=NULL)
## CLV1000 Resin grade visualization


ggplot(Shark_tidy1,aes(y=Global_Power_Avg, x=Ring, color=PRE.DRY)) + geom_point() +coord_cartesian(ylim=
```



##Boxplots (CLV1000) show and seperate the data more than scatter plots.


```
ggplot(Shark_tidy1,aes(y=Global_Power_Avg, x=Ring, fill=PRE.DRY)) + geom_boxplot() +coord_cartesian(yli
```

3.25B SF Jarvis Measurements



```
Shark_tidy2<-filter(Shark.tidy,PRE.DRY=="CLV1000-Yes_Run2"|PRE.DRY=="CLV1000-No_Run2") %>% mutate(Cyl1..
ggplot(Shark_tidy2,aes(y=Global_Power_Avg, x=Ring, fill=PRE.DRY)) + geom_boxplot() +coord_cartesian(ylin
```


3.25B SF Jarvis Measurements



```
## SABIC Resin grade visualization
```

```
Shark_tidy3<-filter(Shark.tidy,PRE.DRY=="0Q1022-Yes_Run2"|PRE.DRY=="0Q1022-No_Run2") %>% mutate(Cyl...0
ggplot(Shark_tidy3,aes(y=Global_Power_Avg, x=Ring, fill=PRE.DRY)) + geom_boxplot() +coord_cartesian(ylin
```


3.25B SF Jarvis Measurements



CLV1000 vs SABIC OQ1022

```
Shark_tidy4<-filter(Shark.tidy,PRE.DRY=="CLV1000-Yes_Run2"|PRE.DRY=="OQ1022-Yes_Run2") %>% mutate(Cyl1..
ggplot(Shark_tidy4,aes(y=Global_Power_Avg, x=Ring, fill=PRE.DRY)) + geom_boxplot() +coord_cartesian(ylin
```

3.25B SF Jarvis Measurements




Facet Wrapping (Rings 1 & 4) to see impact of Power on Target value.

```
Shark.R1<-filter(Shark_tidy1, Ring=="R.01")
Shark.R2<-filter(Shark_tidy1, Ring=="R.02")
Shark.R3<-filter(Shark_tidy1, Ring=="R.03")
Shark.R4<-filter(Shark_tidy1, Ring=="R.04")
Shark.R1.R2.R3.R4<-filter(Shark_tidy1, Ring=="R.01"|Ring=="R.02"|Ring=="R.03"|Ring=="R.04")


ggplot(Shark.R1,aes(y=Global_Power_Avg, x=PRE.DRY, fill=PRE.DRY)) + geom_boxplot() +coord_cartesian(yli
  yintercept = 4.47) + geom_text(aes(label="Target", y=4.46, x=4.4 ), color="red")
```

3.25B SF Jarvis Measurements




```
ggplot(Shark.R4,aes(y=Global_Power_Avg, x=PRE.DRY, fill=PRE.DRY)) + geom_boxplot() +coord_cartesian(ylim=c(4.3, 4.75), yintercept = 4.3) + geom_text(aes(label="Target", y=4.25, x=4.4 ), color="red")
```


3.25B SF Jarvis Measurements



```
ggplot(Shark.R1.R2.R3.R4,aes(y=Global_Power_Avg, x=PRE.DRY, fill=PRE.DRY)) + geom_boxplot() +coord_cart
```

Study 2: 3.25B SF Jarvis Measurements–Drying @120C for 3hrs



Summary:

Inferential statistics: Based on the data visualization and trend analysis, a summary of key variables and observations were pursued. The table below shows a snapshot of describing the statistical summary of all Global Power avg calculations spanning rings 1 through 11. The average value described by the median for each ring was measured to get an idea of its relationship with the theoretical u-lens Power profile. Alongside the median, the standard deviation was derived from the median values to provide an understanding of the spread of the lens data. Furthermore, it is hypothesized that material conditioning has an impact on the optical global power average. The results did consistently show that the condition of PRE DRYING the material had a higher global power average across all 11 rings. In order to see whether this observation is true, **hypothesis testing** using the **students t-test** as well as **ANOVA** testing was utilized. The t-test will compare the difference in means between two samples, while the ANOVA will look at all the different samples in the data (4 total). Test variables included 2 sets of questions.

- Was there a statistically significant difference between global power average for CLV1000 vs SABIC OQ1022 grades?
- Is the difference noticed between DRIED vs NON-DRIED statistically significant using alpha of 0.05?

Hypothesis testing and ANOVA using alpha(0.05) results **did not** show any discrimination between the Avg Power or Radius between the various resin grades, nor did it show that DRYING of the different grades of material made an impact on the optical results. The results were as follows:

- The p-value from the t.test for CLV1000 vs Sabic OQ1022 was 0.9376.(»0.05), Not statistically significant

- The p - value for CLV1000 Dried vs Non-dried was 0.6649. ($\gg 0.05$), Not significant
- F value for ANOVA was 0.105($< F\text{-crit- } 0.9567$) when comparing all different resin and process conditioning, Not significant.


```
library(ggplot2)
library(dplyr)
Summary.shark<-Shark_tidy1 %>%
  group_by(PRE.DRY, Ring) %>%
  summarise(Avg_Power = median(Global_Power_Avg),
            std.dev=sd(Global_Power_Avg)) %>%
  mutate(Radius = round((591-518)/(Avg_Power-(591-518)/(167.81)))) %>%
  mutate(Target = rep(c(4.47,4.52,4.5,4.43,4.3,4.11,3.87,3.58,3.23,2.83,2.37)))
Summary.shark %>%
  filter(PRE.DRY=="CLV1000-No_Run2" | PRE.DRY=="CLV1000-Yes_Run2" | PRE.DRY=="OQ1022-No_Run2" |
         PRE.DRY=="OQ1022-Yes_Run2") %>%
  slice(1:2)
```

```
## # A tibble: 8 x 6
## # Groups:   PRE.DRY [4]
##   PRE.DRY      Ring Avg_Power  std.dev Radius Target
##   <chr>      <chr>    <dbl>    <dbl>  <dbl>  <dbl>
## 1 CLV1000-No_Run2 R.01      4.45    0.0218   18.2   4.47
## 2 CLV1000-No_Run2 R.02      4.47    0.0240   18.1   4.52
## 3 CLV1000-Yes_Run2 R.01      4.61    0.0489   17.5   4.47
## 4 CLV1000-Yes_Run2 R.02      4.61    0.0453   17.5   4.52
## 5 OQ1022-No_Run2  R.01      4.45    1398.    18.2   4.47
## 6 OQ1022-No_Run2  R.02      4.45    0.0390   18.2   4.52
## 7 OQ1022-Yes_Run2 R.01      4.59    0.0570   17.6   4.47
## 8 OQ1022-Yes_Run2 R.02      4.60    0.0573   17.5   4.52
```

```
## Jarvis Clear-reader optical power profile.
```

```
Shark_tidy1 %>%
  group_by(PRE.DRY, Ring) %>%
  summarise(Avg_Power = median(Global_Power_Avg),
            std.dev=sd(Global_Power_Avg)) %>%
  mutate(Radius = round((591-518)/(Avg_Power-(591-518)/(167.81)))) %>%
  mutate(Target = rep(c(4.47,4.52,4.5,4.43,4.3,4.11,3.87,3.58,3.23,2.83,2.37
))) %>%
  ggplot(aes(x=Ring, y=Avg_Power, color=PRE.DRY, group=PRE.DRY)) + geom_point() +geom_line(aes(x=Ring, y=Avg_Power))
  ggtitle("Summarized Data- Jarvis Cleareader Optical Power profile")
```


Summarized Data– Jarvis Cleareader Optical Power profile



Standard deviation vs power profile

```
Shark_tidy1 %>%
  group_by(PRE.DRY, Ring) %>%
  summarise(Avg_Power = median(Global_Power_Avg),
            std.dev=sd(Global_Power_Avg)) %>%
  mutate(Radius = round((591-518)/(Avg_Power-(591-518)/(167.81))) %>%
  mutate(Target = rep(c(4.47,4.52,4.5,4.43,4.3,4.11,3.87,3.58,3.23,2.83,2.37
))) %>%
  ggplot(aes(x=std.dev, y=Avg_Power, color=PRE.DRY, group=PRE.DRY, size=PRE.DRY)) + geom_point() + coord
  ggtitle("Summarized Data- Jarvis Cleareader Optical Power profile")
```

Summarized Data– Jarvis Cleareader Optical Power profile



```
## Facet wrap varied by Ring position
```

```
Shark_tidy1 %>%
```

```
  group_by(PRE.DRY, Ring) %>%
```

```
  summarise(Avg_Power = median(Global_Power_Avg),  
            std.dev=sd(Global_Power_Avg)) %>%
```

```
  mutate(Radius = round((591-518)/(Avg_Power-(591-518)/(167.81)))) %>%
```


```
  mutate(Target = rep(c(4.47,4.52,4.5,4.43,4.3,4.11,3.87,3.58,3.23,2.83,2.37
```

```
))) %>%
```

```
  ggplot(aes(x=std.dev, y=Avg_Power, color=PRE.DRY, group=PRE.DRY, size=PRE.DRY)) + geom_point() + coord
```

```
  ggtitle("Summarized Data– Jarvis Cleareader Optical Power profile") + facet_wrap(~Ring, nrow=3)
```

Summarized Data– Jarvis Cleareader Optical Power profile



```
## Inferential Statistics- t-test
## CLV1000 vs SABIC OQ1022
Sample.data<-Summary.shark %>%
  spread(PRE.DRY, Avg_Power)
t.test(Sample.data$`CLV1000-No_Run2`,Sample.data$`OQ1022-No_Run2`, na.action="na.pass", conf.level = 0.95)

##
## Welch Two Sample t-test
##
## data: Sample.data$`CLV1000-No_Run2` and Sample.data$`OQ1022-No_Run2`
## t = -0.017211, df = 19.997, p-value = 0.9864
## alternative hypothesis: true difference in means is not equal to 0
## 95 percent confidence interval:
## -0.6498991 0.6392627
## sample estimates:
## mean of x mean of y
## 3.783182 3.788500

##Impact of Drying on Power profile
t.test(Sample.data$`CLV1000-Yes_Run2`,Sample.data$`OQ1022-Yes_Run2`, na.action="na.pass", conf.level = 0.95)

##
## Welch Two Sample t-test
##
## data: Sample.data$`CLV1000-Yes_Run2` and Sample.data$`OQ1022-Yes_Run2`
```

```
## t = 0.079281, df = 19.999, p-value = 0.9376
## alternative hypothesis: true difference in means is not equal to 0
## 95 percent confidence interval:
## -0.6477312 0.6989130
## sample estimates:
## mean of x mean of y
## 3.922682 3.897091
```

```
t.test(Sample.data$`CLV1000-Yes_Run2`, Sample.data$`CLV1000-No_Run2`, na.action="na.pass", conf.level = 0.95)
```

```
##
## Welch Two Sample t-test
##
## data: Sample.data$`CLV1000-Yes_Run2` and Sample.data$`CLV1000-No_Run2`
## t = 0.43971, df = 19.967, p-value = 0.6649
## alternative hypothesis: true difference in means is not equal to 0
## 95 percent confidence interval:
## -0.5223486 0.8013486
## sample estimates:
## mean of x mean of y
## 3.922682 3.783182
```

```
##Inferential statistics- ANOVA.
```

```
anova(lm(Summary.shark$Avg_Power~Summary.shark$PRE.DRY))
```

```
## Analysis of Variance Table
##
## Response: Summary.shark$Avg_Power
##              Df Sum Sq Mean Sq F value Pr(>F)
## Summary.shark$PRE.DRY 3  0.173  0.05767   0.105 0.9567
## Residuals          40 21.964  0.54911
```

```
anova(lm(Summary.shark$Radius~Summary.shark$PRE.DRY))
```

```
## Analysis of Variance Table
##
## Response: Summary.shark$Radius
##              Df Sum Sq Mean Sq F value Pr(>F)
## Summary.shark$PRE.DRY 3    7.24    2.414  0.0604 0.9803
## Residuals          40 1598.38   39.960
```

Conclusion & Reflection

The objective of this analysis was to firstly determine whether injection molded spectacle lens for myopia control prevention can be successfully made to conform to certain specification, and secondly, whether different process conditioning(DRY vs NonDRY) had any influence on the optical properties of the spectacle based lens. The graphical results do show that both grade of Polycarbonate resin(**CLV & Sabic**) were within the design specification (+/- 2%), regardless of process condition. Although the data for the CLV1000 grade had a slightly narrower distribution than the Sabic rade when comparing their respective standard deviations, as highlighted by the size change in the **standard deviation vs Power profile** figure above. The narrow distribution is better in terms of process repeatability and reproducibility.

Additionally there does appear to be a Diopter (optical power) difference between DRIED vs NonDied grades of Polycarbonate by roughly (**~0.20D diopters more on avg for DRIED**). The implications of having a slightly higher optical power reading means better vision correction and a higher efficacy for myopia control in children. Interestingly enough, t-test and ANOVA testing did not pick up a the difference in diopter increase since it was so small, and considered it to be non- statistically significant. This is expected from a mathematical standpoint especially when the significance level is so low (0.05), however from a physics and optics standpoint a 0.25D increase in optical power is quite extraordinary. In conclusion, what the results show is that PRE.DRYING of the resin would be the better solution for optimal efficacy when it comes to this technological solution for myopia.

```
knitr::write_bib(x = c("rmarkdown", "knitr"), file = "bib.bib")
```

References

Shao, Zheng. Optical changes during normal emmetropization, lens-induced myopia and its recovery in the young chick eye. MS thesis. University of Waterloo, 2015.

Tarutta, Elena, et al. "Myopia: Why Study the Mechanisms of Myopia? Novel Approaches to Risk Factors Signaling Eye Growth-How Could Basic Biology Be Translated into Clinical Insights? Where Are Genetic and Proteomic Approaches Leading? How Does Visual Function Contribute to and Interact with Ametropia? Does Eye Shape Matter? Why Ametropia at All?." Optometry and vision science 88.3 (2011): 404-447.